



TALENT SURGICAL

PRECISION SURGICAL
INSTRUMENTS CATALOG

TALENT SURGICAL - PREMIUM INSTRUMENT CATALOG

PRECISION SURGICAL TWEEZERS



Precision Surgical Tweezers
Stainless Steel
Serrated Grip
Premium Quality Instrument

ADSON FORCEPS



Adson Tissue & Dressing Forceps
Stainless Steel
1x2 Teeth or Smooth Jaws
Premium Quality Instrument

TALENT SURGICAL

TALENT SURGICAL



Precision Iris Scissors
Stainless Steel Blades
Matte Black Handles
Premium Surgical Instrument

Precision Iris Scissors - TC-100

Designed for delicate procedures, our Premium Iris Scissors feature razor-sharp stainless steel blades and ergonomic matte black handles for superior control and reduced glare. Ideal for ophthalmology, dermatology, and fine dissection. CE Certified. Reusable and Autoclavable.

Blade Length: 11.5 cm

Material: German Stainless Steel & Polymer Coating

Finish: Satin Black Handles, Polished Blades

ADD TO CATALOG REQUEST

TALENT SURGICAL

Premium Surgical Instruments



www.buysurgicalinstruments.com.au | info@buysurgicalinstruments.com.au
Call : +61430855424

TALENT SURGICAL PRECISION INSTRUMENTS



- *Length: 14.5 cm
- *Type: Metzenbaum Style
- *Coating: PVD Black Matt
- *Material: High-Grade Stainless Steel

TALENT SURGICAL

PREMIUM TWO-TONE INSTRUMENT SERIES

Crile-Wood Needle Holders, Hemostats & Suture Scissors







TALENT SURGICAL

Premium Stainless Steel Hemostatic Forceps

Model: TS-SF14 | Straight | Length: 14cm | Fine Serrated Jaws



Model: TS-CF14 | Curved | Length: 14cm | Fully Serrated Jaws

TS

TALENT SURGICAL



Precision Surgical Tweezers
Stainless Steel
Serrated Grip
Premium Quality Instrument

www.buysurgicalinstruments.com.au | info@buysurgicalinstruments.com.au
Call : +61430855424

TALENT SURGICAL



Precision Iris Scissors
Stainless Steel Blades
Gold Plated Handles
Premium Surgical Instrument

www.buysurgicalinstruments.com.au | info@buysurgicalinstruments.com.au

Call : +61430855424

TALENT SURGICAL



Precision Iris Scissors
Stainless Steel Blades
Matte Black Handles
Premium Surgical Instrument

www.buysurgicalinstruments.com.au | info@buysurgicalinstruments.com.au

Call : +61430855424

TALENT SURGICAL



Precision Hemostat Forceps
Stainless Steel
Premium Surgical Instrument

www.buysurgicalinstruments.com.au | info@buysurgicalinstruments.com.au
Call : +61430855424



PREMIUM TITANIUM-COATED STITCH SCISSORS

Model: TS-SS-45-B

Length: 4.5 inches (11.5 cm)

Tips: Precision Sharp/Sharp, curved

Autoclavable: Yes (up to 134°C)

Packaging: Individual Sterile Pouch

TALENT SURGICAL



Gold-Plated Hemostat



Standard Stainless Steel Hemostat



Micro-Line Hemostat



*Gold & Stainless Steel Hemostats
Precision Surgical Scissors
Premium Quality Instruments*